

Product Specification

Product Name: [Digital Thermometer](#) JT002NM

Power Supply: 1.5 DC, 1LR41 Button Battery.

Power Consumption: 0.1MW

Battery Life: About 1100 readings (If used for 10 minutes per day)

Sensing Unit: Thermistor

Measurement Method: Actual measurement

Temperature Display: 3-digit, +°C(°F)display in 0.1 degree increments

Measurement Accuracy: $\pm 0.1^{\circ}\text{C}$ (32°C - 43°C)

$\pm 0.18^{\circ}\text{F}$ (89.60°F - 109.40°F)

Measurement Time: 2m

Measurement Range: 32°C to 43°C (89.6°F to 109.40°F)

Usage Environment: Surrounding temperature: $+10^{\circ}\text{C}$ to $+40^{\circ}\text{C}$ (50°F to 104°F)

Storage Environment: -20°C to $+60^{\circ}\text{C}$ (-4°F to 140°F)

Temp and Humidity: Relative humidity 30-95%RH

Weight: Approx. 15g(battery included)

Dimensions:20mm(W)*130mm(l)*13mm(d)

Accessories: Test battery(Alkaline-Magnesium Button Battery LR4),Storage case,User manual.

Packing: 80pcs/CTN

CTN Size:331mm*305mm*260mm

show details

01 battery cover

This product USES the button battery, push the battery cover and replace the battery



02 Internal Battery

Open the battery cover to see the internal details determine the quality



03 Elliptic switch

A button switch, convenient and concise




04 LED display screen

The LED screen shows the digital value is more intuitive



05 Heat detector

Using soft head probe, safe and convenient



Product measure methods:

axillary measurement

Place the [temperature](#) probe in the middle of the armpit and make the display part of the inner axillary central part of the temperature become high.

Use the arm gently to press the heat probe to fix the sensor under the arm. Push up slightly, clamp, and hold your arm gently with your hand.

Oral method

Stick to the left or right of the base of the tongue

Use your tongue to press the thermometer

Don't open your mouth while trying to avoid the thermometer

Anorectal measurement

Gently insert the probe end of the [thermometer](#) into the rectum for several millimeters. If you feel any discomfort during insertion, stop the insertion immediately

If you want to measure anal temperature for a newborn baby, put your child on a hard bed, lift their legs and maintain a 90-degree Angle to the body so that they don't roll over when measured

A variety of test temperature methods



Oral Thermometry



Oral Thermometry



axillary measurement



Anal measurement method

