

## Non Touch Baby Kopf Thermometer

Eigenschaften:

1. Adopts zuverlässige sensor made in Deutsch.
2. All-neu & patentiertes Sondendesign sorgt für hohe Genauigkeit.
3. Ausgezeichnete Anpassung an Umgebungsbedingungen Temperatur. Präzise und zuverlässig auch in komplizierter Umgebung.
4. Body-Modus und Objekt-Modus verfügbar, ein Schalter zum Modus wechseln
5. Beeper für Fieber oder hohe Temperatur (Alarmwert kann vom Benutzer definiert werden).
6. Automatischer Rückruf der letzten Lesung.
7. Groß LCD mit Hintergrundbeleuchtung.
8. C / F lesen verfügbar.
9. Automatische Abschaltung.

<b>Technical parameters</b>		
Measurement method	Non-contact	
Effective distance	5cm~8cm	
Range	Human body	35.5°C~42.9°C (95.9°F~109.2°F)
	Object	0°C~118°C (32°F~244°F)
Accuracy	Human body	±0.2°C/0.4°F
	Object	±1.0°C/1.8°F
Resolution	0.1°C/°F	
Working condition	10°C~40°C(50°F~104°F) RH≤80%Non-condensing	
Storage condition	-25°C~55°C (-13°F~131°F) RH≤90%Non-condensing	
Power supply	d.c. 3V 2 AA alkaline Battery	
Power consumption	When off≤10uW	
	When measurement≤30mW	
Power level indicator	Indication for low power level	
Memory	Automatic recall of last reading	
Display	LCD back-lit display	
Reading scale	Celsius or Fahrenheit	
Automatic shut off	In 30 seconds	
Dimensions	93mm×153mm×41mm	
Net weight	125g	
Standards	EN60601-1,EN 12470-5 , ASTM 1965-98	

<b>Packing information:</b>	
Packing method:	Each in a color box with plastic tray
Qty/ Ctn:	48pcs
G.W:	13 kg
N.W:	12 kg
Carton Dimension:	56*39*33 CM

<b>Certificates:</b>	ISO13485,CE0499,RoHS
----------------------	----------------------





## Non Contact Infrared Forehead Thermometer




**CE**   **RoHS**   **ISO13485**



Error less than or equal to  $\pm 0.2^{\circ}\text{C}$

1s accurate measure



Kontakt: [sanieleung@dor-medex.com](mailto:sanieleung@dor-medex.com)

Webseite: [Thermometerversorger.com](http://Thermometerversorger.com)