

Medexpro 3d printing technology, please contact us for more information.

3d printing technology, please contact us for more information.



3D printing technology, please contact us for more information.

The filtration effect of non oil particulate matter was at least 96.7%

customizable



Custom items

The custom product itself belongs to the particularity. After confirming the customized sample, the product cannot be changed. It is suggested to confirm the information with the customer service before production.

如何正确佩戴口罩

步骤 1

检查口罩，确认口罩的密封性。将口罩放在手边，不要用手触摸口罩的正面，避免污染。

步骤 2

将口罩的耳带挂在耳朵上。

步骤 3

调整口罩，使其紧密贴合面部。

步骤 4

检查口罩是否完全覆盖口鼻，并确保没有漏气。如果有漏气，请重新调整。

步骤 5

佩戴口罩后，避免用手触摸口罩的正面。如需调整，请使用干净的手指触摸耳带。使用后，将口罩放入密封袋中，避免污染。

