

## ear and forehead baby thermometer

1. 4 in 1: ear measurement, forehead measurement, ambient [temperature](#) and clock display.
2. Back-lit for easy reading.
3. Measurement finished in one second. The perfect integration of advanced [infrared](#) sensing technology and accurate computing ensure it's accuracy and reliability.
4. All functions completed by one button.
5. Memorizes last reading.
6. [Fever alarm](#) by buzzer & icon.
7. With plastic case and storage.

## Ear and Forehead Thermometer



**CE**    **RoHS**    **ISO13485**

### Specification:

Measurement range: 32.0°C-42.9°C (89.6°F-109.2°F)

Allowance:  $\pm 0.2^{\circ}\text{C}/0.4^{\circ}\text{F}$     35.5°C-42.0°C (95.9°F-107.6°F)

$\pm 0.3^{\circ}\text{C}/0.5^{\circ}\text{F}$     32.0°C-35.4°C (89.6°F-95.7°F)

$\pm 0.3^{\circ}\text{C}/0.5^{\circ}\text{F}$     42.1°C-42.9°C (107.8°F-109.2°F)

Resolution: 0.1°C/0.1°F

Operation condition: 10.0°C-40.0°C(50.0°F-104.0°F) Humidity≤80% non-condensing

Storage condition: -25.0°C-55.0°C (-13.0°F-131.0°F) Humidity≤90% non-condensing

Battery: D.C.3V 1×CR2032 Lithium battery

Dimension: 114mm×34mm×26mm

Net weight: Approx.50g (including battery)

Packing info.:

Qty/ctn: 100pcs

Ctn size: 50\*34.5\*33.5cm

G.W./Ctn: 12.00kgs


Remarks:

- 1.Casing colors available: green,pink,blue or customized;
- 2.CE0120 and RoHS compliant.
- 3.Supply capability: 20,000pcs per month.



Error less than or equal to  $\pm 0.2^{\circ}\text{C}$

Is accurate measure



**Contact:** [sanieleung@dor-medex.com](mailto:sanieleung@dor-medex.com)

**Website:** [thermometersupplier.com](http://thermometersupplier.com)