


Unique new design for human Four layers Dust Allergy Flu [Earloop Mask](#) Breath Healthy Safety [Respirator](#) Warm Ski Cycling Half Face Mouth Masks Filters Dust Pollen Allergens Flu Germs Masks



**Features:**

Human body engineering 3D stereo design  
Multilayer non-woven fabric and meltblown, total 4 layers  
Aluminum belt hold shape



### Applicable scope

can be used to protect non oily particles, reduce the odor of some organic compounds.

Non oily particulates: dust (mineral property, metallic dust fiber dust) pollen and haze.

Which can be used to protect environment air pollution fine particle PM 2.5.

Can't be suitable for:

Respirators do not produce oxygen, not in oxygen-deficient environment.

**remember:** Can only be used for pollutants, the wrong use will lead to disease or even death.

professional  
removed dust

## Invisible potential threats in the air



Industrial dust



Sand invasion



Pollen allergens



PM2.5



Secondhand smoke



City abandoned



journey



Harmful bacteria

### Steps for usage

step 1

Open the mask, nose clip out. Both hands were pulled side of the band, so that the nose folder above the mask.

step 2

Put the mask against the chin.

step 3

Pull the band to the ear.

step 4

Place your hands in the middle of the metal nose clip. From the middle of the nose by the status quo, until close to the bridge of the nose.

step 5

Into the work area□Feel the air leak from the nose, should be re-adjust the nose clip.If you feel no air overflow can enter the work area.

