

## **HW-2S [Temassız kızılötesi termometre](#)**

32 MEMORY

TARAMA FONKSİYONU (GERÇEK ZAMANLI [SICAKLIK BODY / OBJECTS MODU](#))

1. insan vücudu için Aralığı: 35.5C-42.9C
2. Nesne için menzil: 0C-118C
3. Arkadan aydınlatmalı ekran
4. C / F mevcut okuma
5. son okuma Otomatik hatırlama
6. Otomatik kapanma
7. Pil seviyesi göstergesi
8. Sapma ayarı: 0C-3.0C
9. Ateş alarmı için buzzer

Uygulanabilir standart: EN 12470-5 ve ASTM 1965-98

Ölçüm mesafesi: 5CM-8CM

Insan vücudu için doğruluk:

$\pm 0.3 \text{ } ^\circ\text{C} / 0.5 \text{ } ^\circ\text{F}$  35.5-42.9  $^\circ\text{C}$  (95.9-107.6  $^\circ\text{F}$ )

Nesne için doğruluk:

$\pm 0.3 \text{ } ^\circ\text{C} / 0.5 \text{ } ^\circ\text{F}$  32.0-35.4  $^\circ\text{C}$  (89.6-95.7  $^\circ\text{F}$ )

Çözünürlük: 0.1  $^\circ\text{C} / 0.1 \text{ } ^\circ\text{F}$

Çalışma koşulu: 10.0  $^\circ\text{C}$ -40.0  $^\circ\text{C}$  (50.0  $^\circ\text{F}$ -104.0  $^\circ\text{F}$ )

Nem% 80 yoğunlaşmayan

Depolama koşulu: -25.0  $^\circ\text{C}$ -55.0  $^\circ\text{C}$  (-13.0  $^\circ\text{F}$ -131.0  $^\circ\text{F}$ )

Nem% 90 yoğuşmasız

Pil: 2 \* AA pil

# FOUR FUNCTIONS IN ONE

Not only measuring body temperature, more care of the whole family's safety



Measuring body temperature



Measuring water temperature



Test milk temperature



Measuring room temperature

Web sitesi: [thermometersupplier.com](http://thermometersupplier.com)

İletişim: [sanieleung@dor-medex.com](mailto:sanieleung@dor-medex.com)